

Common Patterns seen on the OAT

Steve Imgrund, MS, CNS, LDN



Objectives

Big picture overview of OAT profile

Investigate common patterns within the test

Explore other nuanced patterns

What additional testing can be considered?

Organic Acids

WHAT ARE ORGANIC ACIDS (OA)?

Organic acidic products of cellular metabolism that are excreted in urine (in mammals)

- Produced by living organisms including humans, bacteria, and fungi
- Evaluation of these downstream products of metabolic pathways provides insight into potential nutrient deficiencies, inflammation, toxicity, and other imbalances that could be contributing to clinical complaints

Origins of Organic Acid Testing

- To rule out rare Inborn Errors of Metabolism (IEM)- usually in infancy
- Elevations of these organic acids (OA) reflect dysfunction in specific metabolic pathways
 - Accumulation of these toxic metabolites can be life-threatening
 - Symptoms observed in the newborn period include poor feeding and weight gain; nausea and vomiting; neuromuscular issues (e.g., poor tone, seizures); and susceptibility to infection

Use of OAs have evolved from investigating IEM to providing insight into functional metabolic imbalances.

Overview: OAT Fundamentals

Requisition #:

Practitioner:

Patient Name:

Date of Collection:

Date of Birth:

Patient Age:

76

Time of Collection:

Patient Sex:

Print Date:

Report Date:



Organic Acids Test - Nutritional and Metabolic Profile

Metabolic Markers in Urine Reference Range (mmol/mol creatinine) Patient Value Reference Population - Females Age 13 and Over

Intestinal Microbial Overgrowth

Yeast and Fungal Markers

1 Citramalic	≤ 3.6	1.1	
2 5-Hydroxymethyl-2-furoic (Aspergillus)	≤ 14	1.2	
3 3-Oxoglutaric	≤ 0.33	0.09	
4 Furan-2,5-dicarboxylic (Aspergillus)	≤ 16	4.7	
5 Furancarboxylicglycine (Aspergillus)	≤ 1.9	0.17	
6 Tartaric (Aspergillus)	≤ 4.5	0.23	
7 Arabinose	≤ 29	11	
8 Carboxycitric	≤ 29	27	
9 Tricarballic (Fusarium)	≤ 0.44	0.10	

Bacterial Markers

10 Hippuric	≤ 613	35	
11 2-Hydroxyphenylacetic	0.06 - 0.66	0.38	
12 4-Hydroxybenzoic	≤ 1.3	0.46	
13 4-Hydroxyhippuric	0.79 - 17	7.7	
14 DHPPA (Beneficial Bacteria)	≤ 0.38	0.14	

Clostridia Bacterial Markers

15 4-Hydroxyphenylacetic (C. difficile, C. stricklandii, C. lituseburense & others)	≤ 19	15	
16 HPPHA (C. sporogenes, C. caloritolerans, C. botulinum & others)	≤ 208	133	
17 4-Cresol (C. difficile)	≤ 75	1.3	
18 3-Indoleacetic (C. stricklandii, C. lituseburense, C. subterminale & others)	≤ 11	1.8	

Fungus

Non-Specific
Bacteria

Clostridia

Oxalates

Mitochondria Function

Neurotransmitter Pathways

Mosaic Diagnostics

Requisition #:

Practitioner:

Patient Name:

Date of Collection:

Metabolic Markers in Urine Reference Range (mmol/mol creatinine) Patient Value Reference Population - Females Age 13 and Over

Oxalate Metabolites

19	Glyceric	0.77 - 7.0	2.2	
20	Glycolic	16 - 117	75	
21	Oxalic	6.8 - 101	66	

Glycolytic Cycle Metabolites

22	Lactic	≤ 48	14	
23	Pyruvic	≤ 9.1	0.35	

Mitochondrial Markers - Krebs Cycle Metabolites

24	Succinic	≤ 9.3	1.3	
25	Fumaric	≤ 0.94	0.29	
26	Malic	0.06 - 1.8	0.63	
27	2-Oxoglutaric	≤ 35	19	
28	Aconitic	6.8 - 28	7.2	
29	Citric	≤ 507	245	

Mitochondrial Markers - Amino Acid Metabolites

30	3-Methylglutaric	≤ 0.76	0.35	
31	3-Hydroxyglutaric	≤ 6.2	3.3	
32	3-Methylglutaconic	≤ 4.5	1.6	

Neurotransmitter Metabolites

Phenylalanine and Tyrosine Metabolites				
33	Homovanillic (HVA) <i>(dopamine)</i>	0.80 - 3.6	3.4	
34	Vanillylmandelic (VMA) <i>(norepinephrine, epinephrine)</i>	0.46 - 3.7	1.3	
35	HVA / VMA Ratio	0.16 - 1.8	H 2.6	
36	Dihydroxyphenylacetic (DOPAC) <i>(dopamine)</i>	0.08 - 3.5	H 4.4	
37	HVA/ DOPAC Ratio	0.10 - 1.8	0.78	
Tryptophan Metabolites				
38	5-Hydroxyindoleacetic (5-HIAA) <i>(serotonin)</i>	≤ 4.3	1.7	
39	Quinolinic	0.85 - 3.9	3.5	
40	Kynurenic	≤ 2.2	1.1	

Pyrimidines

Ketones

Fatty Acid Oxidation

B Vitamins and Antioxidants

Mosaic Diagnostics

Requisition #:

Practitioner:

Patient Name:

Date of Collection:

Metabolic Markers in Urine Reference Range (mmol/mol creatinine) Patient Value Reference Population - Females Age 13 and Over

Pyrimidine Metabolites - Folate Metabolism

41 Uracil	≤ 9.7	4.3	
42 Thymine	≤ 0.56	0.17	

Ketone and Fatty Acid Oxidation

43 3-Hydroxybutyric	≤ 3.1	0.49	
44 Acetoacetic	≤ 10	0.15	
45 Ethylmalonic	0.44 - 2.8	1.2	
46 Methylsuccinic	0.10 - 2.2	0.87	
47 Adipic	0.04 - 3.8	0.52	
48 Suberic	0.18 - 2.2	1.0	
49 Sebacic	≤ 0.24	0.05	

Nutritional Markers

Vitamin B12			
50 Methylmalonic *	≤ 2.3	0.78	
Vitamin B6			
51 Pyridoxic (B6)	≤ 34	0.54	
Vitamin B5			
52 Pantothenic (B5)	≤ 10	1.1	
Vitamin B2 (Riboflavin)			
53 Glutaric *	0.04 - 0.36	H 0.40	
Vitamin C			
54 Ascorbic	10 - 200	78	
Vitamin Q10 (CoQ10)			
55 3-Hydroxy-3-methylglutaric *	0.17 - 39	6.8	
Glutathione Precursor and Chelating Agent			
56 N-Acetylcysteine (NAC)	≤ 0.28	0	
Biotin (Vitamin H)			
57 Methylcitric *	0.19 - 2.7	1.0	

* A high value for this marker may indicate a deficiency of this vitamin.

Detoxification

Amino Acids

Minerals

Mosaic Diagnostics

Requisition #:

Practitioner:

Patient Name:

Date of Collection:

Metabolic Markers in Urine Reference Range (mmol/mol creatinine) Patient Value Reference Population - Females Age 13 and Over

Indicators of Detoxification

Metabolic Marker	Reference Range (mmol/mol creatinine)	Patient Value	Reference Population - Females Age 13 and Over
Glutathione			
58 Pyroglutamic *	10 - 33	24	
Methylation, Toxic exposure			
59 2-Hydroxybutyric **	0.03 - 1.8	0.81	
Ammonia Excess			
60 Orotic	0.06 - 0.54	0.26	
Aspartame, salicylates, or GI bacteria			
61 2-Hydroxyhippuric	≤ 1.3	0.26	

* A high value for this marker may indicate a Glutathione deficiency.

** High values may indicate methylation defects and/or toxic exposures.

Amino Acid Metabolites

62 2-Hydroxyisovaleric	≤ 2.0	0.07	
63 2-Oxoisovaleric	≤ 2.1	0.06	
64 3-Methyl-2-oxovaleric	≤ 2.0	0.05	
65 2-Hydroxyisocaproic	≤ 2.0	0	
66 2-Oxoisocaproic	≤ 2.0	0.01	
67 2-Oxo-4-methylbutyric	≤ 2.0	0.38	
68 Mandelic	≤ 2.0	0.07	
69 Phenylactic	≤ 2.0	0.02	
70 Phenylpyruvic	≤ 2.0	0	
71 Homogentisic	≤ 2.0	0.03	
72 4-Hydroxyphenylactic	≤ 2.0	0.23	
73 N-Acetylaspartic	≤ 38	1.0	
74 Malonic	≤ 9.7	4.8	
75 4-Hydroxybutyric	≤ 4.8	1.3	

Mineral Metabolism

76 Phosphoric	1,000 - 5,000	2,166	
---------------	---------------	-------	--

The Organic Acids Test (OAT)

Mosaic Diagnostics' OAT provides a comprehensive nutritional and metabolic snapshot of an individual's overall health.

These OAs provide insight into the functional capability of key physiologic areas: [gut health, mitochondrial function, neurotransmitter status, detoxification, and macronutrient breakdown and nutritional status](#)

By understanding these potential contributing factor to several diseases, and by correcting these imbalances, [improvement in symptoms and optimization of the system may be achieved.](#)

This **urine-based assay** assesses 76 organic acids grouped into clinically relevant functional classes:

Intestinal Microbial Overgrowth	Ketone and Fatty Acid Oxidation
Oxalate Metabolites	Nutritional Markers
Glycolytic Cycle Metabolites	Indicators of Detoxification
Mitochondrial Markers / Krebs Cycle Metabolites	Amino Acid Metabolism
Mitochondrial Markers / Amino Acid Metabolites	Mineral Metabolism
Neurotransmitter Metabolites	Indicator of Fluid Intake (Creatinine)
Pyrimidine Metabolites / Folate Metabolism	

OAT Patterns

Pattern recognition – Common findings

Clostridia and Neurotransmitter Metabolism

Yeast activity and Oxalates

Toxic Exposure



Clostridia and Neurotransmitter Metabolism

Clostridia Overgrowth

- Many different clostridia species – Some problematic, but not all
- Clostridia organic acids found on the OAT are known to have potential adverse health effects
 - **4-hydroxyphenylacetic** (*C. difficile*, *C. stricklandii*, *C. lituseburense* & others)
 - **HPHPA** (*C. sporogenes*, *C. caloritolerans*, *C. botulinum* & others)
 - **4-Cresol** (*C. difficile*)
 - **3-Indoleacetic** (*C. stricklandii*, *C. lituseburense*, *C. subterminale* & others)



CDC. (n.d.). Public domain - computer-generated image of a cluster of barrel-shaped, *Clostridium perfringens* bacteria. Retrieved April 2, 2024,


from <https://www.rawpixel.com/image/2288362/free-photo-image-bacterium-bacteria-clostridium-perfringen>

DBH Inhibition


HPHPA

4-cresol


4-HydroxyPhenylacetic



Dopamine







Norepinephrine



Epinephrine (Adrenaline)






Clostridia Bacterial Markers

15	4-Hydroxyphenylacetic <i>(C. difficile, C. stricklandii, C. lituseburense & others)</i>	≤ 30	H 40	
16	HPHPA <i>(C. sporogenes, C. caloritolerans, C. botulinum & others)</i>	≤ 227	214	
17	4-Cresol <i>(C. difficile)</i>	≤ 76	1.1	
18	3-Indoleacetic <i>(C. stricklandii, C. lituseburense, C. subterminale & others)</i>	≤ 11	1.8	

This test was developed, and its performance characteristics determined by Mosaic Diagnostics Laboratory. It has not been cleared or approved by the US Food and Drug Administration.

Neurotransmitter Metabolites

Phenylalanine and Tyrosine Metabolites

33	Homovanillic (HVA) <i>(dopamine)</i>	≤ 14	H 15	
34	Vanillylmandelic (VMA) <i>(norepinephrine, epinephrine)</i>	0.87 - 5.9	3.9	
35	HVA / VMA Ratio	0.12 - 3.0	H 3.9	
36	Dihydroxyphenylacetic (DOPAC) <i>(dopamine)</i>	0.07 - 4.0	H 5.7	
37	HVA / DOPAC Ratio	1.5 - 2.6	H 2.7	

Clostridia Bacterial Markers

15	4-Hydroxyphenylacetic <i>(C. difficile, C. stricklandii, C. lituseburense & others)</i>	≤ 18		8.9	
16	HPHPA <i>(C. sporogenes, C. caloritolerans, C. botulinum & others)</i>	≤ 102	H	159	
17	4-Cresol <i>(C. difficile)</i>	≤ 39		4.9	
18	3-Indoleacetic <i>(C. stricklandii, C. lituseburense, C. subterminale & others)</i>	≤ 6.8		0.30	

Neurotransmitter Metabolites

Phenylalanine and Tyrosine Metabolites

33	Homovanillic (HVA) <i>(dopamine)</i>	0.39 - 2.2	H	3.2	
34	Vanillylmandelic (VMA) <i>(norepinephrine, epinephrine)</i>	0.53 - 2.2		0.75	
35	HVA / VMA Ratio	0.32 - 1.4	H	4.2	
36	Dihydroxyphenylacetic (DOPAC) <i>(dopamine)</i>	0.27 - 1.9	H	2.9	
37	HVA/ DOPAC Ratio	0.17 - 1.6		1.1	

Tryptophan Metabolites

38	5-Hydroxyindoleacetic (5-HIAA) <i>(serotonin)</i>	≤ 2.9		1.3	
39	Quinolinic	0.52 - 2.4		1.0	
40	Kynurenic	≤ 1.8		0.15	

Clostridia Bacterial Markers

15	4-Hydroxyphenylacetic <i>(C. difficile, C. stricklandii, C. lituseburense & others)</i>	≤ 19	17	
16	HPHPA <i>(C. sporogenes, C. caloritolerans, C. botulinum & others)</i>	≤ 208	87	
17	4-Cresol <i>(C. difficile)</i>	≤ 75	3.1	
18	3-Indoleacetic <i>(C. stricklandii, C. lituseburense, C. subterminale & others)</i>	≤ 11	7.5	


This test was developed, and its performance characteristics determined by Mosaic Diagnostics Laboratory. It has not been cleared or approved by the US Food and Drug Administration.

Neurotransmitter Metabolites


Phenylalanine and Tyrosine Metabolites

33	Homovanillic (HVA) <i>(dopamine)</i>	0.80 - 3.6	H 7.7	
34	Vanillylmandelic (VMA) <i>(norepinephrine, epinephrine)</i>	0.46 - 3.7	2.6	
35	HVA / VMA Ratio	0.16 - 1.8	H 2.9	
36	Dihydroxyphenylacetic (DOPAC) <i>(dopamine)</i>	0.08 - 3.5	H 7.3	
37	HVA/ DOPAC Ratio	0.10 - 1.8	1.1	


Yeast and Fungal Markers



Bacterial Markers



Clostridia Bacterial Markers



Neurotransmitter Metabolites

Phenylalanine and Tyrosine Metabolites

33 Homovanillic (HVA) <i>(dopamine)</i>	0.80 - 3.6	H 13	
34 Vanillylmandelic (VMA) <i>(norepinephrine, epinephrine)</i>	0.46 - 3.7	1.3	
35 HVA / VMA Ratio	0.16 - 1.8	H 10.0	
36 Dihydroxyphenylacetic (DOPAC) <i>(dopamine)</i>	0.08 - 3.5	H 6.2	
37 HVA/ DOPAC Ratio	0.10 - 1.8	H 2.1	


Tryptophan metabolites

38 5-Hydroxyindoleacetic (5-HIAA) <i>(serotonin)</i>	≤ 4.3	0.74	
39 Quinolinic	0.85 - 3.9	3.1	
40 Kynurenic	≤ 2.2	0.21	

Neurotransmitter Nuances

- Nuanced Pattern: Non-Clostridia Factors

High HVA/VMA Ratio



PMID:
15060114
23675164
3607034
8515683
988175
24711750
<https://doi.org/10.1136/jmg.2003.009282>

Yeast activity and Oxalates

Production

Candida, and other yeast, can produce excess glyoxylate, liver has the potential to convert that to Oxalate

Molds have been shown to produce Oxalic Acid, directly

Yeast activity and Oxalates

Clinical Overview

Crystal-like structures associated with symptom of pain

- Pain without Etiology?


Three metabolites measured - Oxalic, glycolic, glyceric acids

- Focus on Oxalates


Important to evaluate due to:

- Connection to candida and mold exposure
- Association with B6 need, and/or mitochondrial dysfunction

Glyoxylate Cycle in Yeast




Aspergillus sp (niger) & Oxalate



Intestinal Microbial Overgrowth

Yeast and Fungal Markers




Oxalate Metabolites




Intestinal Microbial Overgrowth

Yeast and Fungal Markers



Oxalate Metabolites




Yeast/Candida Nuances


Intestinal Microbial Overgrowth




Yeast and Fungal Markers









Bacterial Markers



Clostridia Bacterial Markers



Oxalate Metabolites

19	Glyceric	0.77	- 7.0	H 22	 
20	Glycolic	16	- 117	H 128	 
21	Oxalic	6.8	- 101	H 308	 

Nuances for high oxalate markers


- Threshold for Tolerance and Elimination is highly variable
- Address 1st page (if warranted) for major source
- B6 is needed to favor Oxalate to Glycine pathway
- Levels excreted not necessarily reflective of body burden.
- Low Oxalate Diet
 - Dependent on situation, and not dictated by numerical values on test
 - Rarely warranted or mild elevations
 - Evaluate symptom profile for severity



Yeast Activity and Riboflavin

Common pattern linked to acetaldehyde detoxification


- Ethanol is an end product of yeast metabolism of glucose
- Acetaldehyde produced in the previous step
 - Requires NAD, B2, Fe, Molybdenum



Vitamin B2 (Riboflavin)

53 Glutaric *

0.04 - 0.36 H 0.45



The dietary modification and treatment of intestinal *Candida* overgrowth – a pilot study

S Otašević ¹, S Momčilović ², M Petrović ³, O Radulović ⁴, N M Stojanović ³,
V Arsić-Arsenijević ⁵

Affiliations + expand

PMID: 30166063 DOI: [10.1016/j.mycmed.2018.08.002](https://doi.org/10.1016/j.mycmed.2018.08.002)

Abstract

Objective: The aim of this study was to evaluate the effectiveness of an alternative treatment in a form of recommended diet modification during and after conventional treatment with antifungals in

“Results of this pilot study showed that patients who adhered to diet modification during and after treatment with nystatin had better outcomes of ICOG and strongly suggest the need for diet modification in these patients...”

symptomatic effect in 56 out of 80 (70.0%) patients and 29 out of 40 (72.5%) in CG, with no statistically significant difference. However, at the second control stool examination, significantly higher percent (85%) of cured patients was recorded after three months of the recommended diet comparing with CG-17 out of 40 (42.5%).

Conclusion: Results of this pilot study showed that patients who adhered to diet modification during and after treatment with nystatin had better outcomes of ICOG and strongly suggest the need for diet modification in these patients which recommendation could reduce excessive prescription of antifungals.

Keywords: Antifungal treatment; Dietary modification; Intestinal candidosis.



Copyright © 2018 Elsevier Masson SAS. All rights reserved.

[PubMed Disclaimer](#)

Toxic Exposure – Heavy Metals, Mold, Chemicals

TOXIC PATTERNS ON OAT





Glycolytic Cycle Metabolites

22	Lactic	≤ 48		28	
23	Pyruvic	≤ 9.1		3.5	

Mitochondrial Markers - Krebs Cycle Metabolites

24	Succinic	≤ 9.3	H	16	
25	Fumaric	≤ 0.94	H	1.5	
26	Malic	0.06 - 1.8	H	3.6	
27	2-Oxoglutaric	≤ 35		18	
28	Aconitic	6.8 - 28		19	
29	Citric	≤ 507	H	702	

Mitochondrial Markers - Amino Acid Metabolites

30	3-Methylglutaric	≤ 0.76	H	1.3	
31	3-Hydroxyglutaric	≤ 6.2	H	8.8	
32	3-Methylglutaconic	≤ 4.5		3.1	

Succinic acid

Elevations often associated with Fatigue and Brain Fog


Succinate is formed in the Krebs cycle from succinyl coA

In the body succinate is converted to fumarate via succinate dehydrogenase


This enzyme is the only one to participate in the Krebs cycle and ETC

Common concern with elevations: brain fog and fatigue

Look for possible toxins




Mitochondrial Markers - Krebs Cycle Metabolites



Pyroglutamic Acid Production

Gamma-Glutamyl Cycle

- Negative Feed back loop for Glutathione to inhibit gamma- glutamyl Cysteine synthase
- Without glutathione, the feedback inhibition is removed, and there is an upregulation of this enzyme, resulting in increased amounts of gamma-glutamyl cysteine, which itself is a precursor to the organic acid 5-oxoproline (pyroglutamic acid)



Indicators of Detoxification


Glutathione

58 Pyroglutamic * 10 - 33 **H** 42




Methylation, Toxic exposure

59 2-Hydroxybutyric ** 0.03 - 1.8 1.7



Ammonia Excess

60 Orotic 0.06 - 0.54 0.25



Aspartame, salicylates, or GI bacteria


61 2-Hydroxyhippuric ≤ 1.3 **H** 2.1



2-Hydroxybutyric

Homocysteine has 2 choices

- Methionine via Methylation
- Cysteine via trans-sulfuration pathway



1. NCI Thesaurus. (n.d.). Ncithesaurus.nci.nih.gov. Retrieved April 2, 2024, from

https://ncithesaurus.nci.nih.gov/ncitbrowser/ConceptReport.jsp?dictionary=NCI_Thesaurus&ns=ncit&code=C120002

2. Hellman, N. (2009b, June 5). 5-Oxoprolinuria as a Cause for Metabolic Acidosis. Renal Fellow Network. [https://www.renalfellow.org/2009/06/05/5-](https://www.renalfellow.org/2009/06/05/5-oxoprolinuria-as-cause-for-metabolic/)


[oxoprolinuria-as-cause-for-metabolic/](https://www.renalfellow.org/2009/06/05/5-oxoprolinuria-as-cause-for-metabolic/)

Indicators of Detoxification

Glutathione

58 Pyroglutamic *

5.7 - 25 16



Methylation, Toxic exposure

59 2-Hydroxybutyric **


≤ 1.2 H 9.1



Ammonia Excess

60 Orotic


≤ 0.46 0.10



Aspartame, salicylates, or GI bacteria

61 2-Hydroxyhippuric

≤ 0.86 0.19



Case Example

60-year-old female

- Chief complaint - Recurrent yeast infections
- Presenting with joint pain, skin rash and metabolic dysfunction and disrupted sleep
 - High TGs, CRP, hbA1c
- Previously treated with both herbals and nystatin

Intestinal Microbial Overgrowth

Yeast and Fungal Markers

1 Citramalic	≤ 3.6	0.81	
2 5-Hydroxymethyl-2-furoic <i>(Aspergillus)</i>	≤ 14	1.5	
3 3-Oxoglutaric	≤ 0.33	H 0.67	
4 Furan-2,5-dicarboxylic <i>(Aspergillus)</i>	≤ 16	1.3	
5 Furancarboxylglycine <i>(Aspergillus)</i>	≤ 1.9	0.64	
6 Tartaric <i>(Aspergillus)</i>	≤ 4.5	H 7.0	
7 Arabinose	≤ 29	H 334	
8 Carboxycitric	≤ 29	0.17	
9 Tricarballic <i>(Fusarium)</i>	≤ 0.44	0.12	
Bacterial Markers			
10 Hippuric	≤ 613	144	
11 2-Hydroxyphenylacetic	0.06 - 0.66	0.44	
12 4-Hydroxybenzoic	≤ 1.3	H 6.9	
13 4-Hydroxyhippuric	0.79 - 17	3.4	
14 DHPPA (Beneficial Bacteria)	≤ 0.38	0.13	
Clostridia Bacterial Markers			
15 4-Hydroxyphenylacetic <i>(C. difficile, C. stricklandii, C. lituseburense & others)</i>	≤ 19	13	
16 HPPHA <i>(C. sporogenes, C. caloritolerans, C. botulinum & others)</i>	≤ 208	2.9	
17 4-Cresol <i>(C. difficile)</i>	≤ 75	37	
18 3-Indoleacetic <i>(C. stricklandii, C. lituseburense, C. subterminale & others)</i>	≤ 11	7.9	

Oxalate Metabolites

19	Glyceric	0.77 - 7.0	H 8.5	
20	Glycolic	16 - 117	48	
21	Oxalic	6.8 - 101	H 243	

Glycolytic Cycle Metabolites

22	Lactic	≤ 48	12	
23	Pyruvic	≤ 9.1	2.4	

Mitochondrial Markers - Krebs Cycle Metabolites

24	Succinic	≤ 9.3	H 9.4	
25	Fumaric	≤ 0.94	H 1.2	
26	Malic	0.06 - 1.8	H 2.3	
27	2-Oxoglutaric	≤ 35	19	
28	Aconitic	6.8 - 28	15	
29	Citric	≤ 507	507	

Mitochondrial Markers - Amino Acid Metabolites

30	3-Methylglutaric	≤ 0.76	0.32	
31	3-Hydroxyglutaric	≤ 6.2	4.8	
32	3-Methylglutaconic	≤ 4.5	1.7	

Neurotransmitter Metabolites

Phenylalanine and Tyrosine Metabolites

33	Homovanillic (HVA) <i>(dopamine)</i>	0.80 - 3.6	H 8.7	
34	Vanillylmandelic (VMA) <i>(norepinephrine, epinephrine)</i>	0.46 - 3.7	1.3	
35	HVA / VMA Ratio	0.16 - 1.8	H 6.6	
36	Dihydroxyphenylacetic (DOPAC) <i>(dopamine)</i>	0.08 - 3.5	H 8.1	
37	HVA / DOPAC Ratio	0.10 - 1.8	1.1	

Pyrimidine Metabolites - Folate Metabolism

41 Uracil	≤ 9.7	2.9	
42 Thymine	≤ 0.56	0.13	

Ketone and Fatty Acid Oxidation

43 3-Hydroxybutyric	≤ 3.1	H 3.8	
44 Acetoacetic	≤ 10	5.2	
45 Ethylmalonic	0.44 - 2.8	1.2	
46 Methylsuccinic	0.10 - 2.2	2.0	
47 Adipic	0.04 - 3.8	2.8	
48 Suberic	0.18 - 2.2	H 6.3	
49 Sebacic	≤ 0.24	H 0.63	

Nutritional Markers


Vitamin B12			
50 Methylmalonic *	≤ 2.3	1.4	
Vitamin B6			
51 Pyridoxic (B6)	≤ 34	3.4	
Vitamin B5			
52 Pantothenic (B5)	≤ 10	8.8	
Vitamin B2 (Riboflavin)			
53 Glutamic *	0.04 - 0.36	H 0.96	
Vitamin C			
54 Ascorbic	10 - 200	L 2.5	
Vitamin Q10 (CoQ10)			
55 3-Hydroxy-3-methylglutaric *	0.17 - 39	26	
Glutathione Precursor and Chelating Agent			
56 N-Acetylcysteine (NAC)	≤ 0.28	0.02	
Biotin (Vitamin H)			
57 Methylcitric *	0.19 - 2.7	1.0	

Indicators of Detoxification

Glutathione

58 Pyroglutamic *


10 - 33 21



Methylation, Toxic exposure

59 2-Hydroxybutyric **


0.03 - 1.8 H 2.1



Ammonia Excess

60 Orotic

0.06 - 0.54 0.39



Aspartame, salicylates, or GI bacteria

61 2-Hydroxyhippuric

≤ 1.3 0.10



Pattern recognition – Common findings

Clostridia and Neurotransmitter Metabolism

Yeast activity and Oxalates

Toxic Exposure

A wide-angle photograph of a long, straight asphalt road stretching towards the horizon. The road is flanked by a dry, brown desert landscape. The sky is a deep blue with wispy white clouds, and a bright sun is shining in the upper left, creating a lens flare effect. The overall mood is one of vastness and forward movement.

An OAT test is just the start



Complementary Testing with Other Mosaic Diagnostics Profiles

Glyphosate Test

- Clostridia, mitochondrial stress

MycoTOX

- Fungal activity, mitochondrial stress, detoxification

Metals – Toxic + Nutrient Elements (Hair vs Urine vs Blood)

- Mitochondrial stress, detoxification section, NT section

IgG FoodMap * Candida

- Microbial section

Complementary Testing with Other Mosaic Diagnostics Profiles

Comprehensive Stool Analysis

- Microbial section - Bacterial analytes

Saliva Hormone Profile

- NT section

DNA Methylation Pathway Profile

- NT section, detoxification section

Thank You

Q & A